## Taber®Rotary Abraser 5135/5155

#### Wear and Abrasion Testing

Taber abrading wheels produce characteristic rub wear action.

Mounted to a rotating turntable, specimens are subjected to the rub-wear action of two abrasive wheels. Driven by the test sample, the wheels produce abrasion marks that form a pattern of crossed arcs over a circular ring approximately 30 square centimeters. This reveals abrasion resistance at all angles relative to the weave or grain of the material.

### Referenced in numerous international test standards and specifications.

Simple to operate, this instrument has been accepted worldwide as a standard for wear and abrasion research, quality and process control, materials evaluation, and product development. Taber Rotary Abrasers are capable of providing reliable data in a matter of minutes, compared to the years that may be required with in-use testing.

#### Test a wide range of products.

The Taber Abraser can be used to test virtually any flat specimen. Its field of application has included tests of: painted, lacquered, powder coated, and electroplated surfaces; textile fabrics ranging from sheer silks to heavy upholstery and carpeting; and solid materials such as metals, stone and ceramics. Other materials include paper, glass, plastics, leather, rubber, linoleum, laminates, plus many others.

5155 ABRASER

Model 5155 Shown

AS9100B & ISO 9001:2008 CE

# Taber® Rotary Abraser 5135/5155

Wear and Abrasion Testing

